

ActivPal Use Guidelines

Instructions for Application

1. Pick a leg that you would like the ActivPal to be on. Swab the skin in the center of your thigh with an alcohol wipe and allow it to dry.
2. Peel off the “Hypafix” adhesive layer on the backside of the ActivPal.
3. Stick the wrapped ActivPal on your leg. The arrow should be facing **UP**.
4. Check daily that the device is secure and in place. If the ActivPal becomes loose, or if skin irritation occurs, see the following instructions on how to re-wrap and apply.



Re-Applying Your ActivPal

1. Carefully cut through layers of tape to remove the rubber-wrapped device.
2. Place the ActivPal in the center of your Hypafix square on the fabric side. Arrow should be facing **UP**.
3. Take one piece of green tape. Carefully remove layer labeled “Remove First” and adhere directly on the device and the Hypafix, across the black arrow.



Stay close to the surface of the device and the Hypafix, there should be **NO** air gaps under the tape. Once stuck, peel off green film.

4. Adhere a second piece of green tape along the arrow, forming a cross. Follow directions from step 3. Remove green film when finished.



5. See Instructions for Application above.